

Technical Data Sheet

10.09

Product:	PUDOL Safe & Clean
Article number:	8190
Product type:	Moist disinfecting wipes
Composition:	Cleaning wipes with disinfecting solution.
pH-Value:	pH-Value = ab. 7
Characteristics:	For disinfection of surfaces in all areas, where an immediate disinfection is required. PUDOL Safe & Clean features safe disinfection properties, pleasant fragrance and a good material compatibility. Dries quickly and clean.
Application:	Open the pail, pull out the damp cloth and tear it at the perforation. Close the cover again. Suitable for medical and dental clinics, hospitals, kindergartens, schools, public areas and food establishments, wherever a rapid disinfection is required, f. ex. treatment chairs, baby's changing tables, bedsteads, furniture and general furniture's such as treatments couches, cabinets, appliances, doors and sanitary facilities.; f. ex. toilet seats, bathtubs, etc. Reaction time: 1 minute Don't apply on surfaces which are not sensitive against alcohol, f. ex. acrylic glass.
Effect:	Bactericidal, fungicidal, tuberculocidal and limited virucidal (incl. HIV, HBV, HCV, Influenza-viruses like H1N1, H5N1) (acc. DVV and RKI)
Ingredients (acc. EEC-recommend.):	Alcohols, quaternary ammonium compounds and auxiliaries.
Packaging:	70 cloths / pail 4 or 8 pail / crt.

This information is without any obligation. We reserve conditional modifications by product improvement or necessary raw material exchange. Users must check the product suitability on not mentioned materials on their own authority. Therefore, our advices are without any obligation and liability claims cannot be assert against us.